



LABFACILITY


TEMPERATURE & PROCESS TECHNOLOGY

Tipo T Cable de Termopar Rizado RetrÃfÃ¡ctil (ANSI)



XE-7034-001

TYPE-T (ANSI)



XE-7034-001

TYPE-T (ANSI)



XE-7034-001

TYPE-T (ANSI)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Disponible en tipo T

Color codificado en ANSI

2 metros de longitud cuando se extienden por completo

Terminado con un enchufe en miniatura instalado

Exterior de poliuretano

Especificaciones

Specifications

Product Code	XE-7034-001
Cable Length	2M fully extended
Cable Type	Polyurethane outer
Connector Type	Thermocouple Plug
Connector Size	Miniature
Colour	Blue
Standards Met	Colour code (cores & jacket) to ANSI MC96.1
Thermocouple Type	T
Termination Type	Miniature Thermocouple Plug