



# LABFACILITY


TEMPERATURE & PROCESS TECHNOLOGY

## Tipo T Cable de Termopar Rizado RetrÃfÃctil (ANSI)



XE-7034-001

### TYPE-T (ANSI)



XE-7034-001

### TYPE-T (ANSI)



XE-7034-001

### TYPE-T (ANSI)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Disponible en tipo T

Color codificado en ANSI

2 metros de longitud cuando se extienden por completo

Terminado con un enchufe en miniatura instalado

Exterior de poliuretano

### **Especificaciones**

## **Specifications**

<b>Product Code</b>	XE-7034-001
<b>Cable Length</b>	2M fully extended
<b>Cable Type</b>	Polyurethane outer
<b>Connector Type</b>	Thermocouple Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Blue
<b>Standards Met</b>	Colour code (cores & jacket) to ANSI MC96.1
<b>Thermocouple Type</b>	T
<b>Termination Type</b>	Miniature Thermocouple Plug