




6 núcleos aislados con PFA, pantalla de cobre chapada en estaño, cable / cable del sensor PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

6 núcleos aislados con PFA, pantalla de cobre chapada en estano, cable / cable del sensor PRT

Alta temperatura a 260°C, teflón PFA extruido® construcción. Comúnmente utilizado para termómetros Pt100?, también adecuado con otros tipos como Pt130, Pt500 y Pt1000?

Características técnicas**Specifications**

General Description	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores with nylon cord, tin plated copper braided screen with PFA insulated outer jacket
Cable Length	10, 25m, 50m
Cable Type	PFA Teflon® construction
Core / Strands	7/0.2mm
Insulation Material	PFA Insulated
Max. Temperature	+260°C
Min. Temperature	-75°C
Number of Wires	6
Standards Met	Cores colour coded Red & White in accordance with IEC-751, black outer jacket