



**6 núcleos aislados con PFA, pantalla de cobre chapada en estaño, cable / cable del sensor PRT**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**6 núcleos aislados con PFA, pantalla de cobre chapada en estano, cable / cable del sensor PRT**

Alta temperatura a 260°C, teflón PFA extruido® construcción. Comúnmente utilizado para termómetros Pt100?, también adecuado con otros tipos como Pt130, Pt500 y Pt1000?

**Características técnicas****Specifications**

<b>General Description</b>	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores with nylon cord, tin plated copper braided screen with PFA insulated outer jacket
<b>Cable Length</b>	10, 25m, 50m
<b>Cable Type</b>	PFA Teflon® construction
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	PFA Insulated
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Number of Wires</b>	6
<b>Standards Met</b>	Cores colour coded Red & White in accordance with IEC-751, black outer jacket