




# 4 núcleos aislado de PVC, pantalla de cobre chapada en estaño, cable / cable del sensor PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**4 núcleos aislado de PVC, pantalla de cobre chapada en estano, cable / cable del sensor PRT**

Cable aislado de PVC para extender termómetros de resistencia de platino. Comúnmente utilizado para termómetros Pt100? pero también adecuado con otros tipos como Pt130, Pt500 y Pt1000?

**Características técnicas****Specifications**

<b>General Description</b>	4 core PVC insulated 7/0.2mm tin plated copper conductors, twisted cores, polyester tape & nylon cord, tin plated copper braided screen with PVC insulated outer jacket
<b>Cable Length</b>	10, 25, 50m
<b>Cable Type</b>	Mylar shield with tin plated copper braiding
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	PVC insulation
<b>Max. Temperature</b>	+105°C
<b>Min. Temperature</b>	-10°C
<b>Number of Wires</b>	4
<b>Standards Met</b>	Cores colour coded 2 x Red / 2 x White to IEC-751 with black outer jacket