

## Accesorios de compresin de acero inoxidable - rosca cnica (BSPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Para adaptarse a una amplia gama de sensores con aislamiento mineral y otros tipos de sonda ofrecidos por Labfacility.

Los accesorios de compresin de acero inoxidable con roscas BSPT (British Standard Pipe Taper) se utilizan ampliamente para conectar tuberías y tubos en diversos entornos exigentes.

Las roscas BSPT son cónicas, lo que significa que el diámetro de la rosca disminuye hacia el final, lo que mejora el sellado a medida que se aprieta el accesorio.

El uso de acero inoxidable proporciona una excelente resistencia a la corrosión, durabilidad y resistencia, lo que hace que estos accesorios sean ideales para las duras condiciones que se encuentran en industrias como el procesamiento químico, el petróleo y el gas, y las aplicaciones marinas.

### Características técnicas

## Specifications

<b>General Description</b>	Male in-line stainless steel compression type with nut & olive
<b>Material</b>	Stainless steel, BSPT tapered thread
<b>Thread</b>	1/8"BSPT, 1/4"BSPT, 1/2"BSPT, 3/8"BSPT or 3/4"BSPT available
<b>Probe Diameter</b>	Suitable for both metric and imperial probe diameters from 1.0mm to 12.0mm