

Accesorios de compresin de latn - rosca cnica (BSPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Para adaptarse a una amplia gama de sensores con aislamiento mineral y otros tipos de sonda ofrecidos por Labfacility.

Los accesorios de compresin de latn con roscas BSPT (British Standard Pipe Taper) se utilizan ampliamente para conectar tuberias y tubos en diversas aplicaciones.

Las roscas BSPT son cnicas, lo que significa que el dimetro de la rosca disminuye hacia el final, mejorando el sellado a medida que se aprieta el accesorio.

Fabricados en latn, estos accesorios proporcionan una buena resistencia a la corrosin, durabilidad y facilidad de instalacin, lo que los hace adecuados para plomera, HVAC y sistemas de presin baja a media que involucran agua, aire y gas.

Caractersticas tcnicas

Specifications

General Description	Male in-line stainless steel compression type with nut & olive
Material	Brass
Thread	1/8"BSPT, 1/4"BSPT, 3/8"BSPT or 1/2"BSPT available
Probe Diameter	Suitable for both metric and imperial probe diameters from 1.5mm