

## Reductores - Latn



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Los reductores de latn son componentes vitales en los sistemas de tuberias, ya que facilitan la transicin entre diferentes tamaos de tubera para cumplir con los requisitos especficos de la aplicacin.

Su resistencia a la corrosin, durabilidad y fuerza los hacen adecuados para una amplia gama de entornos y aplicaciones.

La seleccin, instalacin y mantenimiento adecuados de los reductores de acero inoxidable garantizan un rendimiento eficiente y confiable.

### Caractersticas tcnicas

## Specifications

<b>General Description</b>	A reducer allows for a change in size to meet specific requirements of the application
<b>Material</b>	Brass
<b>Sizes</b>	BSPP