




## Termopar de unin expuesta de alambre desnudo de calibre fino (conductores de 0,025 mm) - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopar ultrafino de punta soldada tipo K para una sensibilidad ptima.

### Caracterstias tcnicas

## Specifications

<b>Product Code</b>	XE-3517-001
<b>General Description</b>	Very fast response to a fraction of a second, Ideal for minimal invasive temperature sensing.
<b>Sensor Type</b>	Exposed Junction
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	200mm
<b>Cable Type</b>	Fine gauge bare wire
<b>Core / Strands</b>	0.001" (0.025mm)
<b>Max. Temperature</b>	315°C