




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Almohadillas de fijacin trmica en tela de vidrio PTFE



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Almohadillas autoadhesivas en tela de vidrio PTFE para localizar termopares de punta soldada

Clasificacin de -73 C a 260 C

Tamao de la almohadilla 13x18mm

La hoja contiene 20 almohadillas

## Specifications

Product Code	XE-3599-001
--------------	-------------