

Parche autoadhesivo PT100 sensor - tipo PRT - 5M de longitud



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caracterstias tcnicas**Specifications**

Product Code	XE-3661-001
General Description	Suitable for general purpose industrial & commercial use
Sensor Type	Silicone rubber patch with self adhesive foil backing
Element	Pt100
Patch dimensions	L40 x W15mm
Cable Length	5M
Cable Type	PTFE insulated twisted leads
Core / Strands	7/0.2mm
Cable Termination	Bare tails
Termination Type	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
Max. Temperature	150°C
Min. Temperature	-50°C
Accuracy	Class B
Number of Wires	4
Series	Patch Pt100