



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Sonda de penetracin de mano de acero
inoxidable - Tipo K ANSI 6,0 mm x 120 mm**




Image
Coming
Soon

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonda de penetracin de acero inoxidable 316 con punta puntiaguda para la medicin de temperatura de lquidos y semislidos. Compuesto por un asa, un cable en espiral de 2 metros y un mini enchufe.

Necesita calibrar su sensor de temperatura?

Trazable a los estndares de calibracin ISO17025 (UKAS), Labfacility ofrece la opcin de una calibracin trazable de 3 o 5 puntos utilizando nuestra instalacin de calibracin interna, puede seleccionar cualquier punto de temperatura entre -10 C y +200 C.

Tenga en cuenta que si agrega una opcin de calibracin a su pedido, nos pondremos en contacto con usted con respecto a la informacin de entrega.

Caractersticas tcnicas

Specifications

Product Code	XE-3869-001
General Description	Penetration Probe with a pointed tip
Sensor Type	K ANSI
Probe Diameter	6.0mm
Probe Length	120mm
Cable Length	2m
Cable Termination	Miniature Plug
Max. Temperature	+400°C