



## Termopar de mano cargado por resorte con punta de disco de cobre - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Funda de acero inoxidable con mango moldeado y 2 metros de cable de extensin enrollado terminado en un tapn de termopar en miniatura.

## Specifications

<b>Product Code</b>	XE-3816-001
<b>General Description</b>	Surface temperature spring loaded thermocouple with copper disc tip
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Cable Length</b>	2m Fully Extended
<b>Cable Type</b>	Retractable Curly Thermocouple Lead - Polyurethane outer
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	+600°C
<b>Dimensions</b>	Probe 4.7mmO x 63mm long stem, O8 x 17mm tip with O4.5mm copper disc
<b>Tip</b>	Copper disc tip
<b>Series</b>	Handheld Probe