


Tipo K IEC Superficie en ngulo recto



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Funda de acero inoxidable con mango moldeado y 2 metros de cable de extensin enrollado terminado en un tapn de termopar en miniatura.

Specifications

Product Code	XE-3819-001
General Description	Surface temperature fast response thermocouple
Sensor Type	Right angled Probe
Thermocouple Type	K
Calibration	IEC
Cable Length	2m
Cable Termination	Fitted Plug
Max. Temperature	+900°C
Dimensions	Probe O12mm Tip x 40mm Hot Leg and O6mm x 180mm Cold Leg. Coiled element diameter 7mm
Tip	Ceramic Tip, Coiled Element