





# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Termopar de disco de alambre fino IEC, PTFE - Tipo K

XE-3515-001

### FINE WIRE DISC THERMOCOUPLE TYPE K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopar de disco de cobre con tiras adhesivas de kapton de alta temperatura (13x18mm x 5) para una fijacin rpida y fcil a la superficie/proceso.

## Specifications

<b>General Description</b>	Fast response – ideal for surface measurement
<b>Sensor Type</b>	Copper disc, 6.35mm diameter x 0.3mm thick
<b>Thermocouple Type</b>	K
<b>Cable Length</b>	1m
<b>Cable Type</b>	PTFE insulated twisted wire
<b>Core / Strands</b>	1/0.2mm
<b>Termination Type</b>	Bare tails
<b>Max. Temperature</b>	+200°C
<b>Min. Temperature</b>	-0°C
<b>Standards Met</b>	IEC-584-3