


## Termopar de unin expuesta ANSI con aislamiento de fibra de vidrio con enchufe en miniatura - Tipos K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Necesita calibrar su sensor de temperatura?**

Trazable a los estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de una calibración trazable de 3 o 5 puntos utilizando nuestras instalaciones de calibración internas.

**Tenga en cuenta que si agrega una opción de calibración a su pedido, nos pondremos en contacto con usted con respecto a la información de entrega.**

**Características técnicas****Specifications**

<b>General Description</b>	Fast response, welded exposed junction
<b>Sensor Type</b>	Good temperature resistance (not suitable in fluids)
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	300mm, 1, 2, 5 or 10 metre
<b>Cable Type</b>	glassfibre insulated flat pair cable
<b>Core / Strands</b>	1/0.3mm conductors
<b>Cable Termination</b>	Miniature Plug
<b>Max. Temperature</b>	+350°C (short periods up to 400°C)
<b>Min. Temperature</b>	-60°C
<b>Tolerance</b>	ANSI MC96.1 Class 2
<b>Standards Met</b>	ANSI MC96.1