



Sensor de parche autoadhesivo PT100 - Tipo PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Suitable for general purpose industrial & commercial use
Sensor Type	Silicone rubber patch with self adhesive foil backing
Element	Pt100
Patch dimensions	L40 x W15mm
Probe Diameter	6.0mm
Probe Length	Various
Sheath Material	316 stainless steel sheath.
Cable Length	2 or 5m
Cable Type	PTFE insulated twisted leads
Core / Strands	7/0.2mm
Cable Termination	Bare Tails
Termination Type	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
Max. Temperature	+150°C
Min. Temperature	-50°C
Accuracy	class B
Number of Wires	4