



Parche autoadesivo PT1000 Sensor - Tipo PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

General Description	Pt1000 element encapsulated in a silicone rubber patch, Self-adhesive, aluminium foil backing.
Sensor Type	Pt1000 (1000? @ 0°C), thin film Class B to IEC751
Element	Pt1000 sensor to IEC751 Class B
Patch dimensions	Length 30mm x Width 15mm x Height 4mm
Cable Length	1m
Cable Type	Teflon ® insulated twin twisted
Core / Strands	7/0.2mm
Cable Termination	Bare Tails
Max. Temperature	+150°C
Min. Temperature	-50°C