


## Detectores de heridas de alambre de 100 ohmios



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Tenga en cuenta: Los plazos de entrega DIN de clase 1/10 variarán - por favor pregunte o se le contactará con un plazo de entrega después de realizar su pedido.**

### Especificaciones

## Specifications

<b>General Description</b>	Wirewound construction for optimum performance and stability
<b>Sensor Type</b>	Pt100 Single Element
<b>Max. Temperature</b>	+650°C
<b>Min. Temperature</b>	-200°C
<b>Tolerance</b>	Class A, Class B or 1/10 DIN. 1/10DIN relative to IEC 60751 over the range 0°C to 150°C (1/10 of Class B)
<b>Detector Dimensions</b>	Various width / length combination available