

## Calentadores de boquilla Mica



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Estos productos están disponibles hasta agotar existencias con cantidades limitadas disponibles para cada versión, se le avisará al agregar a la cesta si no tenemos suficiente stock - [por favor contáctenos](#) si tiene requisitos específicos.

## Specifications

<b>General Description</b>	Best quality Mica for high dielectric strength and maximum heat transfer
<b>Element</b>	High performance precision wound element for superior heat distribution
<b>Cable Length</b>	1 metre long stainless steel braided cable
<b>Max. Temperature</b>	Recommended for temperatures up to 300°C
<b>Output</b>	240V
<b>Features #1</b>	Full width low expansion straps
<b>Features #2</b>	Specially treated finish to resist oxidation and enhance thermal conductivity
<b>Important Information</b>	Wide range of applications including plastic process machinery, injection / blow moulding, rubber moulding, test equipment, food process industry etc