


## Termopar de autoclave IEC - Dplex tipo T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Esta sonda de carga del esterilizador ofrece un sensor dplex tipo T y es adecuada para su uso en condiciones adversas, incluidos la mayora de los esterilizadores de vapor.

## Specifications

<b>Product Code</b>	XE-9510-001
<b>General Description</b>	First 2.5m section in a flexible stainless steel conduit
<b>Thermocouple Type</b>	Type T Duplex
<b>Thermocouple Junction</b>	2 x insulated tip junctions
<b>Probe Diameter</b>	4mm
<b>Probe Length</b>	250mm
<b>Sheath Material</b>	Stainless steel (316) sheath
<b>Cable Type</b>	Cable is silicon rubber throughout
<b>Max. Temperature</b>	+200°C
<b>Min. Temperature</b>	-50°C