


Termopar de autoclave IEC - Dplex tipo T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Esta sonda de carga del esterilizador ofrece un sensor dplex tipo T y es adecuada para su uso en condiciones adversas, incluidos la mayora de los esterilizadores de vapor.

Specifications

Product Code	XE-9510-001
General Description	First 2.5m section in a flexible stainless steel conduit
Thermocouple Type	Type T Duplex
Thermocouple Junction	2 x insulated tip junctions
Probe Diameter	4mm
Probe Length	250mm
Sheath Material	Stainless steel (316) sheath
Cable Type	Cable is silicon rubber throughout
Max. Temperature	+200°C
Min. Temperature	-50°C