

Caja de equipo ABS negro



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caja ABS negro viene con conjunto de espuma completa que se puede arrancar en cubos de 10 mm y modificado a los requisitos de los usuarios.

Dimensiones

Interno (lxxwxh)
420 x 298 x 131mm

Externo (lxxanto)
432 x 363 x 138mm

Peso (Kg)
1.4

Specifications

Product Code	XE-6050-001
---------------------	-------------