


Insertos Pt100 con funda de cermica - 2 hilos con longitud de cable de 50 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Características técnicas**Specifications**

Product Code	XE-3666-001
General Description	Pt100 thin film element housed in ceramic sheath
Element	Pt100
Cable Length	50mm
Cable Type	PTFE insulated wires
Core / Strands	7/0.15mm
Cable Termination	Bare Tails
Max. Temperature	250°C
Min. Temperature	-50°C
Tolerance	Class B
Detector Dimensions	Ceramic 5mm dia x 35mm long
Number of Wires	2
Series	Pt100 Inserts