



## Insertos Pt100 con funda de ceramica - 4 hilos con longitud de paso de 450 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Características técnicas

**Specifications**

<b>Product Code</b>	XE-3667-001
<b>General Description</b>	Pt100 thin film element housed in ceramic sheath
<b>Element</b>	Pt100
<b>Cable Length</b>	450mm
<b>Cable Type</b>	PTFE insulated wires
<b>Core / Strands</b>	7/0.15mm
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Tolerance</b>	Class B
<b>Detector Dimensions</b>	Ceramic 5mm dia x 35mm long
<b>Number of Wires</b>	4
<b>Series</b>	Pt100 Inserts