


## Termmetro de resistencia higienico Pt100, accesorio estilo RJT de 1,5 " con transmisor incorporado - 3 hilos, 6,0 mm x 125 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Las uniones RJT (tipo de junta de anillo) a menudo se denominan uniones estándar. El compromiso de una junta tipo junta trica que deja una pequeña grieta donde el revestimiento y la parte macho de la unin se superponen.

Ideal para procesos higienicos, como aplicaciones lcteas, cerveceras y farmaceuticas.

### Caractersticas tcnicas

## Specifications

<b>Product Code</b>	XE-5526-001
<b>General Description</b>	Ideal for hygienic processes including dairy, brewery and pharmaceutical applications.
<b>Sensor Type</b>	Pt100, 3 wire single element class A
<b>Head Type</b>	KNEF
<b>Connection Type</b>	Pt100, 3 wire single element class A
<b>Transmitter</b>	SEM-203p transmitter, ranged -50°C to +150°C
<b>Support Dia (D2)</b>	8.0mm
<b>Support Length (L2)</b>	50mm