


**Termometro de resistencia higienico Pt100,  
accesorio Tri-Clamp de 1,5" con transmisor  
integrado - 3 hilos, 6,0 mm x 75 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ideal para procesos higienicos, como aplicaciones lcteas, cerveceras y farmaceuticas.

**Caractersticas tcnicas****Specifications**

|                        |  |
|------------------------|--|
| <b>Product Code</b>    | XE-5529-001                                  |
| <b>Sensor Type</b>     | Pt100, 3 wire single element class A         |
| <b>Probe Diameter</b>  | 6.0mm  |
| <b>Probe Length</b>    | 75mm   |
| <b>Head Type</b>       | KNEF   |
| <b>Transmitter</b>     | SEM-203p transmitter, ranged -50°C to +150°C |
| <b>Number of Wires</b> | 3  |