

**Termometro de resistencia higienico Pt100,
accesorio Tri-Clamp de 1,5" con transmisor
integrado - 3 hilos, 6,0 mm x 75 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ideal para procesos higienicos, como aplicaciones lcteas, cerveceras y farmaceuticas.

Caracterstias tcnicas**Specifications**

Product Code	XE-5529-001
Sensor Type	Pt100, 3 wire single element class A
Probe Diameter	6.0mm
Probe Length	75mm
Head Type	KNEF
Transmitter	SEM-203p transmitter, ranged -50°C to +150°C
Number of Wires	3