




Pt100 4 cable clase B Termómetro de resistencia, cabezal compacto KNS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Process Connection	1/2"BSP x 6mm
Probe Diameter	6mm
Probe Length	Various probe lengths
Sheath Material	316 stainless steel
Head Type	KNS
Max. Temperature	+450°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67
Accuracy	Class B
Number of Wires	4