



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Pt100 Termometro de resistencia de 4 hilos clase B, cabezal KNS compacto con accesorio de 1/2" BSP x 6 mm con longitud de sonda de 100 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Caracterstias tcnicas**Specifications**

Product Code	XE-5509-001
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Process Connection	1/2"BSP x 6mm
Probe Diameter	6.0mm
Probe Length	100mm
Sheath Material	316 Stainless Steel
Head Type	KNS
Max. Temperature	+450°C (connection head @ 170°C)
Min. Temperature	-50°C
IP Rating	IP67
Accuracy	Class B
Number of Wires	4
Series	Pt100 KNS