



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Pt100 Termometro de resistencia de clase B de 4 hilos, cabezal KNS compacto con accesorio de 1/2 "BSP x 6 mm con longitud de sonda de 200 mm**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Caracterstias tcnicas****Specifications**

<b>Product Code</b>	XE-5511-001
<b>Sensor Type</b>	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
<b>Process Connection</b>	1/2"BSP x 6mm
<b>Probe Diameter</b>	6.0mm
<b>Probe Length</b>	200mm
<b>Sheath Material</b>	316 Stainless Steel
<b>Head Type</b>	KNS
<b>Max. Temperature</b>	+450°C (connection head @ 170°C)
<b>Min. Temperature</b>	-50°C
<b>IP Rating</b>	IP67
<b>Accuracy</b>	Class B
<b>Number of Wires</b>	4