

Termmetro infrarrojo IR RayTemp 2



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Uso industrial: precisión no apta para uso médico.

Certificado de calibración trazable gratuito

Registra las temperaturas máximas/mínimas

Muestra la temperatura ambiente

Termómetro infrarrojo RayTemp 2 - Diseado y fabricado aquí en el Reino Unido, el nuevo termómetro infrarrojo sin contacto RayTemp 2. Basta con apuntar el termómetro al objetivo y pulsar el botón de medición para mostrar la temperatura de la superficie en el rango de -49,9 a 349,9 °C con una resolución de 0,1 °C o 1 °C y una precisión garantizada de 1 °C en el rango de 0 a 100 °C.

Certificado de calibración trazable GRATUITO incluido.

Características técnicas

Specifications

Product Code	XE-3941-001
General Description	The RayTemp® 2 incorporates a 5:1 optic ratio (target distance/diameter ratio) and a fixed emissivity of 0.95 making it suitable for a wide range of food and industrial applications.
Type	infrared
Accuracy	±1 °C (0 to 100 °C) otherwise ±2 °C or ±2 % of reading, whichever is greater
Resolution	0.1°C & 1°C
Emissivity	0.95 fixed
Dimensions	25 x 56 x 128mm
Input / Ranges	-49.9 to 349.9°C
Weight	130 grams
Power Consumption	5000 hours continuous use
Power Supply	3 x 1.5 volt AAA



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Important Information	Not suitable for medical use
------------------------------	------------------------------