

## Accesorios de compresin de acero inoxidable - rosca NPT (NPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Para adaptarse a una amplia gama de sensores con aislamiento mineral y otros tipos de sonda ofrecidos por Labfacility.

Los accesorios de compresin de acero inoxidable con roscas NPT (National Pipe Thread) se utilizan ampliamente para conectar de forma segura tuberías y tubos en diversas aplicaciones industriales.

Las roscas NPT son cónicas, lo que significa que el diámetro de la rosca disminuye hacia el final, creando un sello más hermético a medida que se aprieta el accesorio.

El material de acero inoxidable ofrece una excelente resistencia a la corrosión, durabilidad y resistencia, lo que hace que estos accesorios sean adecuados para entornos hostiles, incluido el procesamiento químico, el petróleo y el gas, y las aplicaciones marinas.

### Características técnicas

## Specifications

<b>General Description</b>	Male in-line stainless steel compression type with nut & olive
<b>Material</b>	NPT - Stainless Steel Compression Fitting
<b>Thread</b>	1/8" NPT, 1/4" NPT, 1/2" NPT, 3/4"NPT or 3/8"NPT available
<b>Probe Diameter</b>	Suitable for both metric & imperial probe diameters from 1.0mm to 12.0mm