

Termmetro industrial Therma 1 (Tipo K)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

TENGA EN CUENTA: El plazo de entrega ser de 7-10 das hbiles al elegir la calibracin

El Therma 1 cuenta con pantallas LCD grandes y fciles de leer con circuito abierto 'Err' e indicacin de batera baja. Cada termmetro funciona con tres bateras AAA de 1,5 voltios que dan un mnimo de cinco aos de duracin de la batera. Las unidades se apagan automticamente despues de diez minutos, maximizando la duracin de la batera. Esta funcin puede ser deshabilitada por el usuario, si es necesario.

No se incluye sonda, elija entre la amplia gama de sondas Labfacility para adaptarse a este producto

Los sensores de temperatura porttiles opcionales para usar con este dispositivo estn disponibles en el men desplegable: 26 sondas disponibles, seleccione el estilo, la longitud y el dimetro de la sonda Tipo K.

Caractersticas tcnicas

Specifications

Product Code	XE-3988-001
General Description	Housed in a robust ABS case that contains 'Biomaster' additive that helps to reduce bacterial growth.
Type	Industrial Digital Thermometer
Thermocouple Type	Type K
Accuracy	$\pm 0.4^{\circ}\text{C} \pm 0.1\%$
Input / Ranges	$0.1^{\circ}\text{C} -99.9$ to 299.9°C . 1°C 300 to 1372°C
Resolution of data display	0.1°C & 1°C
Dimensions	25 x 56 x 128mm
Power Supply	Battery 3 x 1.5 volt AAA. Battery life 10000 hours/min 5 years
Sensor Type	Sensor type K miniature thermocouple
Weight	130 grams