


Termopar de pernos IEC, acero inoxidable de fibra de vidrio sobre cable trenzado - Tipo K, J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Los sensores de pernos generalmente se atornillan en la superficie, lo que normalmente aumenta la velocidad de respuesta y la sensibilidad de la medicación. Apto para su uso hasta +250C.

Specifications

Thermocouple Type	K or J
Thermocouple Junction	Bolts have grounded junctions
Cable Length	2m
Cable Type	Glassfibre stainless steel over braided cable
Cable Termination	Bare tails
Max. Temperature	+250°C
Bolt thread/pitch	M8 x 1mm / M10 x 1mm / M12 x 1mm