





LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Termopar de unin expuesta de calibre fino ANSI con aislamiento de PFA (conductores de 0,076 mm) - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Necesita calibrar su sensor de temperatura?

Trazable a los estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de una calibración trazable de 3 o 5 puntos utilizando nuestras instalaciones de calibración internas.

Tenga en cuenta que si agrega una opción de calibración a su pedido, se comunicará con usted con respecto a la información de entrega.

Características técnicas**Specifications**

Product Code	XE-2006-001
General Description	Fine Gauge Exposed Welded Tip Thermocouple
Sensor Type	Fast response, welded tip exposed junction
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	1m
Cable Type	Teflon® PFA insulated wires
Core / Strands	0.003" (0.076mm) conductors
Cable Termination	Bare tails
Max. Temperature	+260°C
Min. Temperature	-75°C
Standards Met	ANSI colour coded