





# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## Termopar de unin expuesta de calibre fino ANSI con aislamiento de PFA (conductores de 0,076 mm) - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Necesita calibrar su sensor de temperatura?**

Trazable a los estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de una calibración trazable de 3 o 5 puntos utilizando nuestras instalaciones de calibración internas.

**Tenga en cuenta que si agrega una opción de calibración a su pedido, nos pondremos en contacto con usted con respecto a la información de entrega.**

**Características técnicas****Specifications**

<b>General Description</b>	Fine Gauge Exposed Welded Tip Thermocouple
<b>Sensor Type</b>	Fast response, welded tip exposed junction
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	300mm, 1m
<b>Cable Type</b>	Teflon® PFA insulated wires
<b>Core / Strands</b>	0.003" (0.076mm) conductors
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Standards Met</b>	ANSI colour coded