

Accesorios de compresin de latn - rosca NPT (NPT)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Para adaptarse a una amplia gama de sensores con aislamiento mineral y otros tipos de sonda ofrecidos por Labfacility.

Los accesorios de compresin de latn con roscas NPT (National Pipe Thread) se usan comunmente para conectar tuberias y tubos en diversas aplicaciones.

Las roscas NPT son cnicas, lo que crea un sello ms hermtico a medida que se aprieta el accesorio.

El latn, conocido por su resistencia a la corrosin, durabilidad y maquinabilidad, hace que estos accesorios sean adecuados para aplicaciones de plomera, HVAC y presin baja a media que involucran agua, aire y gas.

Caractersticas tcnicas

Specifications

General Description	Male in-line stainless steel compression type with nut & olive
Material	Brass
Thread	1/8" or 1/4" NPT available
Probe Diameter	To suite 3.0mm or 6.0mm probe diameters