



# 4 núcleo de goma de silicona aislado, sección redondeada, para sensor cable / cable



25M



**4 CORE**  
XE-2798-025

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**4 núcleos de goma de silicona aislado, sección redondeada, prT sensor cable / cable**

Se puede utilizar para extender hasta 4 configuraciones RTD de alambre. Comúnmente utilizado con termómetros Pt100, pero adecuado con otros tipos de Pt130/500/1000, etc.

**Especificaciones****Specifications**

<b>General Description</b>	4 core silicone rubber insulated 7/0.2mm tin plated copper conductors, twisted together with silicone rubber outer jacket
<b>Cable Length</b>	10, 25m
<b>Cable Type</b>	Silicone rubber insulated
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	Silicone rubber insulated, tin plated copper conductors, twisted together with silicone rubber outer jacket
<b>Max. Temperature</b>	+200°C
<b>Min. Temperature</b>	-40°C
<b>Number of Wires</b>	4
<b>Standards Met</b>	Cores colour coded 2 x Red / 2 x White to IEC-751