

Termmetro infrarrojo IR GM300E (solo para uso no mdico)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Este termómetro infrarrojo se utiliza para medir la temperatura de la superficie del objeto, que es aplicable para varios objetos calientes, peligrosos o difíciles de tocar sin contacto de forma segura y rápida.

Para cantidades mayores, utilice la pestaña "consulta de producto" para ponerse en contacto con nuestro equipo de ventas para obtener una cotización e información sobre el plazo de entrega.

Tenga en cuenta:

estos son termómetros IR industriales, así que verifique las especificaciones antes de comprar, ya que se aplicará una tarifa de reposición del 20% a los productos no utilizados. La precisión de estos artículos los hace inadecuados para uso médico

Características técnicas

Specifications

Product Code	XE-4000-001
Model	GM300E
Accuracy	+/- 1.5°C / +/- 1.5%
Resolution	0.1°C or 0.1°F
Emissivity	0.10-1.00 Adjustable, 0.95 Preset (see download for emissivity chart)
Distance to Spot Size	12:1
Battery	9V
Function #1	MAX / MIN / AVG / DIF reading
Function #2	°C / °F unit selectable
Function #3	Auto Power off & Data Hold
Function #4	High / Low Temperature Alarm setup
Function #5	Laser ON / OFF selectable
Function #6	Backlight ON / OFF selectable

Unit selectable	°C / °F
Power Supply	Requires one standard 9V battery (included)
Measuring range	-50°C to 420°C (-58-788°F)
Important Information	These are industrial IR thermometers, so please check the specifications before purchasing as a restocking fee of 20% will apply to unused products