

GM1311 Termometro digital profesional para alimentos de alto rendimiento



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

TENGA EN CUENTA: El tiempo de entrega ser de 7 a 10 das hbiles al elegir la calibracin

GM1311 Termmetro digital profesional de alto rendimiento para alimentos/barbacoa/carne

Temperatura ambiente :0 hasta 40 C (32 a 104 F)

Humedad ambiental: 10-90% RH

Voltaje de tasa: 3V

1 x Termmetro para alimentos

1 x Manual de usuario en ingls

1 x CR2032

Caractersticas tcnicas

Specifications

Product Code	XE-4002-001
Accuracy	-10 - 100°C: +/-1°C / 100 - 200°C: +/-1.5°C
Resolution	0.1°C / °F
Colour	White
Dimensions	9.29 in x 1.34 in x 0.71 in (23.6 cm x 3.4 cm x 1.8 cm). Display dimension:44x18mm
Function #1	LCD display
Function #2	Temperature lock
Function #3	Max/Min measuring
Function #4	Conversion between Celsius degree and Fahrenheit degree
Function #5	Auto power off
Max. Temperature	300°C (572°F)
Min. Temperature	-50°C (-58°F)
Probe Length	149.5 x 6mm



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Weight	Weight: 0.95 oz (27 g)
Power Supply	CR2032 button battery x 1 (included)