


## Termopar de velcro IEC, PFA con enchufe en miniatura - Tipo K,T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Los termopares de velcro se pueden utilizar como un mtodo rpido y sencillo para medir la temperatura alrededor de la superficie de una tubera.

## Specifications

<b>Thermocouple Type</b>	K or T
<b>Cable Length</b>	1.5M (0.5M Inside Velcro)
<b>Cable Type</b>	PFA
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Miniature flat pin plug, colour coded in accordance IEC 584
<b>Max. Temperature</b>	+150°C
<b>Length</b>	Velcro: 500mm