

## IEC Tipo K 0.7kg Botn de traccin Magnet Termopar, PTFE aislado con enchufe miniatura o estndar



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>General Description</b>	Ideal for ferrous metal surface measurement
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded at tip
<b>Magnet</b>	12.5mm diameter button type magnet. Max pull: 0.7 KG
<b>Cable Length</b>	1 or 2 metres
<b>Cable Type</b>	PTFE insulated twin-twist wires
<b>Cable Termination</b>	Miniature or Standard Fitted Plug
<b>Max. Temperature</b>	+250°C (plug 200°C)