



Termopar de banda magnética IEC tipo K, PTFE aislado con tapn miniatura



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un termopar de banda magnética flexible capaz de unirse a la mayoría de las superficies ferromagnéticas, la gran resistencia del imán garantiza una buena adhesión que permite la simplicidad de la instalación. Está en una amplia gama de aplicaciones de monitoreo de temperatura, particularmente en pruebas y desarrollo.

Specifications

General Description	Ideal for ferrous metal surface measurement
Thermocouple Junction	Grounded measuring junction
Magnet	50 x 25 x 1mm magnetic strip
Cable Length	1 or 2 metre
Cable Type	PTFE insulated wires, twin twisted
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	100°C
Min. Temperature	-50°C