



## Termopar de banda magnética IEC tipo K, PTFE aislado con tapn miniatura



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Un termopar de banda magnética flexible capaz de unirse a la mayoría de las superficies ferromagnéticas, la gran resistencia del imán garantiza una buena adhesión que permite la simplicidad de la instalación. Está en una amplia gama de aplicaciones de monitoreo de temperatura, particularmente en pruebas y desarrollo.

## Specifications

<b>General Description</b>	Ideal for ferrous metal surface measurement
<b>Thermocouple Junction</b>	Grounded measuring junction
<b>Magnet</b>	50 x 25 x 1mm magnetic strip
<b>Cable Length</b>	1 or 2 metre
<b>Cable Type</b>	PTFE insulated wires, twin twisted
<b>Termination Type</b>	Miniature flat pin plug, colour coded in accordance IEC 584
<b>Max. Temperature</b>	100°C
<b>Min. Temperature</b>	-50°C