



Termopar de disco de alambre fino IEC, PTFE con enchufe en miniatura - Tipo K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Tiras adhesivas kapton de alta temperatura (13x18mm x 5) para una fijación rápida y fácil a la superficie/proceso. Disco de cobre de 6,35 mm de diámetro x 0,3 mm de grosor.

Specifications

General Description	Exposed junction thermocouple with Copper Disc brazed at tip
Sensor Type	Copper disc
Thermocouple Type	Type 'K' or 'J' IEC, Type 'K' BS
Thermocouple Junction	Grounded measuring junction
Cable Length	1 or 2 metre
Cable Type	1/0.2mm PTFE insulated twisted wire
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584