


## Termopar de unin expuesta ANSI con aislamiento de PFA con tapn moldeado - Tipo K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Necesita calibrar su sensor de temperatura?**

Trazable a los estándares de calibración ISO17025 (UKAS), Labfacility ofrece la opción de una calibración trazable de 3 o 5 puntos utilizando nuestras instalaciones de calibración internas.

**Tenga en cuenta que si agrega una opción de calibración a su pedido, se comunicará con usted con respecto a la información de entrega.**

**Características técnicas****Specifications**

<b>General Description</b>	Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids, gasses, steam etc.
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Exposed
<b>Cable Length</b>	1, 2 or 3m
<b>Cable Type</b>	PFA Teflon® twin twisted pair
<b>Core / Strands</b>	1/0.3mm
<b>Termination Type</b>	moulded-on type K ANSI miniature thermocouple plug
<b>Max. Temperature</b>	260°C
<b>Min. Temperature</b>	-75°C
<b>Tolerance</b>	Tolerance Class 1 for optimum accuracy
<b>Standards Met</b>	ANSI MC96.1 Colour Coded