



Termopar de bayoneta ajustable IEC, cable trenzado de acero inoxidable de fibra de vidrio - tipo K,J



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Termopar de bayoneta ajustable tipo K o J IEC y tipo J BS

Termopar de bayoneta ajustable con resorte de 200 mm con tapa de fijacin ajustable (11,5 mm de dimetro interno de una sola ranura).

Caracterstias tcnicas**Specifications**

General Description	Thermocouple with adjustable insertion depth
Sensor Type	6.0mm diameter tip
Thermocouple Type	K, J IEC J BS
Thermocouple Junction	Grounded tip for fast response
Cable Length	2, 3m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Cable Termination	Bare wire or standard plug
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	IEC-584-3 / BS
Bayonet cap	11.5mm I.D. single slot
Connector Type	Plug
Connector Size	Standard
Dimensions	25x13mm shim (0.2mm thick)



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thread	M6 x 1mm thread
Spring	200mm