

## Master Class Thermmetro de cocina plegable



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**TENGA EN CUENTA: El plazo de entrega ser de 7-10 das hbiles al elegir la calibracin**

Perfecto para cocinar con precisin carne de res, ternera, cordero, cerdo, pollo y pavo, asegrese de tener un delicioso asado de carne para su almuerzo del domingo, impresionando tanto a amigos como a familiares.

La pantalla LCD fcil de leer muestra claramente la temperatura interna medida por la sonda de acero inoxidable, con un rango de temperatura de -45 C a 200 C y -49 F a 392 F.

Este termmetro plegable tambn es ideal para almacenar fcilmente en los cajones de la cocina, ya que el mecanismo de plegado reduce el tamao total.

### Caractersticas tcnicas

## Specifications

<b>Product Code</b>	XE-3966-001
<b>General Description</b>	The easy to read LCD display clearly shows the internal temperature measured by the stainless steel probe
<b>Sensor Type</b>	Folding Stainless Steel Probe
<b>Max. Temperature</b>	200°C (392°F)
<b>Min. Temperature</b>	-45°C (-49°F)
<b>Dimensions</b>	25.8 x 8.8 x 2.2 cm
<b>Power Supply</b>	1 x LR44 battery (included)
<b>Weight</b>	100g