

GM816 Anemmetro Anemmetro Anemmetro / Anemmetro Digital



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

GM816 Anemmetro Anemmetro Anemmetro / Anemmetro Digital

Ideal para fines industriales y domsticos, como medir la velocidad del viento y la temperatura de la CPU, ventiladores de computadora, aires acondicionados, sopladores de aire, etc.

Specifications

Product Code	XE-4003-001
Model	GM816 Anemometer Wind Meter Anemometer / Digital Anemometer
Accuracy	+/- 5%
Air Temp: Accuracy	+/- 2°C, +/- 3.6F
Air Temp: Range	-10~45°C, 14~113°F
Air Temp: Resolution	Resolution 0.2°C, 0.36°F
Air Velocity: Range	0.3 - 30m/s, 0 - 90km/h, 0 - 5860ft/min, 0 - 65mph, 0 - 55Knots
Air Velocity: Resolution	0.1m/s, 0.3km/h, 19ft/min, 0.2mph, 0.2Knots
Air Velocity: Threshold	0.1m/s, 0.3km/h, 39ft/min, 0.2mph, 0.1Knots
Battery	1 x 3V CR2032 battery (Not Included)
Colour	Yellow / Black
Current Consumption	3mA Approx
Dimensions	10.5 x 4.5 x 2 cm / (4.13 x 1.77 x 0.79)
Operating Humidity	less than or equal to 90%RH
Operating Temperature	-10~45°C (14~113°F)
Weight	105g / 3.7oz
Storage Temperature	-40~60°C, -40~140°F
Thermometer	NTC Thermometer